

Fast help is double help!

Now you don't have to decide which eye to save!



# **Duo Eye Wash**


- Mobile first aid solution for both eyes at the same time.



# **Duo Eye Wash Solutions**

Several eye accidents involve both eyes. Dust, dirt, chemicals and other harmful substances can enter the eye tissue within seconds and in worse case lead to serious eye injuries.

The new Duo Eye Wash Solutions sets new standards for first aid eye wash bottles by offering the ability to rinse both eyes at the same time – without losing valuable time! Duo is therefore a necessary complement to your current eye wash system. You optimize first aid and creates an even higher standard of safety in your workplace. Duo rinses with a gentle even flow while rinsing eye wash fluid and the damaging substances away from the eyes. Furthermore the registered design of the eye device on top of the bottle is supporting the ability to keep the eyes open during rinsing.









#### 1000 ml Plum Eye Wash Duo

Sterile 0.9 % Sodium chloride solution. For injuries with foreign bodies like dust, dirt etc. Complete with an eye device for rinsing both eyes at the same time. Also for refilling of the Duo Eye Wash stations.

Durability: 3 years Rinsing time: approx. 5 minutes Art.-No.: 4800 Units per box: 6



#### 500 ml pH Neutral Duo

Sterile 4.9 % Phosphate buffer solution. For neutralising acids and alkalis . Complete with an eye device for washing both eyes at the same time. Also for refilling of the Duo eye wash station.

**Durability:** 3 years **Rinsing time:** approx. 2 minutes

Art.-No.: 4801 Units per box: 6



#### Plum Eye Wash Station Duo

Wall mounted station containing 1 x 1000 ml Plum Eye Wash Duo bottle, wall mounting device, pictogram and mirror. Especially for general eye washing at smaller work places.

Art.-No.: 4802 Units per box: 1 Station



#### Plum CombiStation Duo

Wall mounted station containing 500 ml pH Neutral Duo and 1000 ml Plum Eye Wash Duo, wall mounting device, pictogram and mirror. Especially suitable in places where there is a risk of eye injuries involving acids and alkalis.

Art.-No.: 4803 Units per box: 1 Station



## Plum Eye Wash box Duo

Dust-proof wall station containing 2 x 1000 ml Plum Eye Wash Duo, pictogram on the front side and a mirror inside the lid. Especially suitable for dirty and dusty production areas and mobile work places.

Art.-No.: 4816 Units per box: 1 Station



### Plum CombiBox Duo

Dust-proof wall station containing 1 x 500 ml pH Neutral Duo and 1 x 1.0 l Plum Eye Wash Duo, pictogram on the front side and a mirror inside the lid. Especially suitable for dusty production areas where there is a risk of eye injuries involving acids and alkali.

Art.-No.: 4810 Units per box: 1 Station

## Fast and simple!









**IMPORTANT:** Rinse with pH Neutral until the bottle is empty. Continue rinsing with Plum Eye Wash until you reach a doctor.



Plum A/S • DK-5610 Assens Tel. +45 6471 2112 • Fax +45 6471 2125 info@plum.dk • www.plum.eu