

# First Aid for institutions, offices, the manufacturing industry, the food industry and the health sector

Make the path to first aid easy, quick and effective. Plums product range includes a first aid assortment developed as complete first aid solutions for any work environment and need. Make the workplace safe and create a sense of security with a customized first aid solution for your company.

## First aid solutions



### QuickFix Elastic

- Elastic and breathable plasters
- Hygienically packed plasters
- Easy to refill and mount

#### Content

Plaster dispenser with 90 Elastic plasters

#### Art. No.

5502

#### Dimensions

W:23 cm x H:13.5 cm



### QuickFix Detectable

- Blue detectable plasters
- Hygienically packed plasters which are particularly breathable
- Easy to refill and mount

#### Content

Plaster dispenser with 90 Detectable plasters

#### Art. No.

5503

#### Dimensions

W:23 cm x H:13.5 cm



### Quick Fix&Rinse

- Quick access to first aid
- Hygienic - individually packed plasters/eye wash in ampoules
- Easy to refill and mount

#### Content

Dispenser with 45 plasters & 5 x 20 ml eye wash ampoules

#### Art. No.

4631

#### Dimensions

W:23 cm x H:13.5 cm



### Plum Eye Wash Station with 2 bottles

- Dust-proof wall station with Plum Eye Wash for use in case of injuries involving foreign bodies, dust and dirt
- Suited for areas with dust and dirt as well as mobile workplaces
- With pictograms and mirror

#### Content

Wall station with 2 x 500 ml Plum Eye Wash

#### Art. No.

4650

#### Dimensions

H:27 cm x W:22.5 cm x D:11 cm

4653 w/heating-device



### Plum Eye Wash Station with pH Neutral and Plum Eye Wash

- Dust-proof wall station with Plum Eye Wash for use in case of injuries involving foreign bodies, dust, dirt, etc. and pH Neutral for injuries involving acids and alkali
- Suited for areas containing a risk for accidents with acids and alkali
- With pictograms and mirror

#### Content

Wall station in polystyrene with 200 ml pH Neutral and 500 ml Plum Eye Wash

#### Art. No.

4787

#### Dimensions

H:27 cm x W:22.5 cm x D:11 cm



## Plum Combibox DUO

- Dust-proof wall station with Plum Eye Wash for injuries involving foreign bodies, dust, dirt, etc. and pH Neutral for acids and alkali
- Suited for areas containing a risk for eye injuries involving acids and alkali in both eyes
- With pictograms and mirror

Content	Art. No.	Dimensions
Wall station in polystyrene with 500 ml pH Neutral DUO and 1000 ml Plum Eye Wash DUO	4810	H:35,7 cm x W:22.5 cm x D:11 cm



## QuickSafe Complete

- Complete and hygienic first aid station containing eye wash, plaster, wound cleanser, wound dressing, burn gel, etc.
- First aid close at hand and for every situation
- Easy to refill and mount

Content	Art. No.	Dimensions
18 QuickCool, 1 pck QuickStop, 20 QuickClean, 2 x 500 ml Eye Wash, 1 x 200 ml pH Neutral, 5 x 20 ml QuickRinse, 45 QuickFix Elastic	5174	H:43 cm x W:25.3 cm x D:9,2 cm

## Refills



### QuickFix Elastic

- Elastic and breathable plasters
- Individually packed

Content	Art. No.	Packaging
45 Elastic plasters	5512	6 x 45 pcs.



### QuickFix Detectable

- Blue detectable and breathable plasters
- Individually packed

Content	Art. No.	Packaging
45 Detectable plasters	5513	6 x 45 pcs.



### QuickRinse refill Eye wash ampoules

- Handy and hygienic ampoules
- Quick first aid for removing foreign bodies, dust and dirt in the eyes

Content	Art. No.	Packaging
5 x 20 ml ampoules (0,9% sterile sodium chloride)	5160	1 x 5 pcs.



### 500 ml Plum Eye Wash Rinsing time: 5 min

- Eye wash for removing foreign bodies, dust and dirt
- Usable for the refill of stations or alone

Content	Art. No.	Packaging
500 ml Plum Eye Wash (0,9% sterile sodium chloride)	4604	12 x 500 ml



### 200 ml pH Neutral Rinsing time: 2 min

- Neutralizes acids and alkali
- Usable for the belt bag, the tool box or similar places as well as for the refill of stations

Content	Art. No.	Packaging
200 ml pH Neutral (4,9% sterile phosphate buffer)	4752	10 x 200 ml



### 500 ml pH Neutral DUO Rinsing time: 2 min

- Neutralizes acid and alkali in both eyes at the same time
- Usable for the refill of stations or alone

Content	Art. No.	Packaging
500 ml pH Neutral (4,9% sterile phosphate buffer)	4801	6 x 500 ml



### 1000 ml Plum Eye Wash DUO Rinsing time: 5 min

- For eye injuries involving foreign bodies, dust and dirt in both eyes
- Usable for the refill of stations or alone

Content	Art. No.	Packaging
1000 ml Plum Eye Wash (0,9% sterile sodium chloride)	4800	6 x 1000 ml